



Item #: EC3600 (MI390)

extriCARE® 3600

Negative Pressure Wound Therapy System

Distributor Sell Sheet



extriCARE® 3600

Negative Pressure Wound Therapy System

The ExtriCARE® 3600 NPWT System is designed to optimize wound healing by effectively removing wound exudates, infectious material, and tissue debris from the wound bed. It administers continuous or intermittent negative pressure to suit individual patient needs. This powerful yet user-friendly device is ideal for various healthcare settings.

Enhancing Wound Care with ExtriCARE®

The ExtriCARE® 3600 utilizes advanced Hybrid hNPWT technology, allowing the use of both traditional foam dressing kits and anatomically fitted bandage dressings. This versatility promotes effective wound healing by providing tailored solutions for various wound types.

Features and Benefits:

- ✓ **Easy to Use:** Our intuitive design and user-friendly interface simplify operation, making it accessible for healthcare professionals and caregivers alike.
- ✓ **Versatile:** Equipped with a high-capacity, high-flow pump, the ExtriCARE® 3600 ensures comprehensive care for a wide range of wound types.
- ✓ **Powerful:** The long-lasting rechargeable battery supports up to 24 hours of continuous treatment, providing reliable performance when needed most.

Order Information

Order Code	Product Descriptions	Packaging
EC3600 (MI390)	extriCARE 3600 NPWT unit	4 units per shipper
EC3600-400C-E	400cc Canister for extriCARE 3600 unit	20 pieces per shipper
EC3600-1000C-E	1000cc Canister for extriCARE 3600 unit	20 pieces per shipper
EC3600-CA400	Carrying case for EC3600 unit utilizing 400cc canister	50 pieces per shipper
EC3600-CA1000	Carrying case for EC3600 unit utilizing 1000cc canister	50 pieces per shipper
EC3600-CHGR	Charger for extriCARE 3600 unit	20 pieces per shipper
EC3600-CLAMP	Bed/Pole clamp for extriCARE 3600 unit	20 pieces per shipper

Device Specifications

Product Code:	EC3600 (MI390)
Weight:	2.87 lbs (1.3kg)
Battery Type:	Lithium (rechargeable)
Vacuum Modes:	Continuous or intermittent
Pressure Range:	40-200 mmHg
Charge Time:	< 5 hours
Canister size:	400cc and 1000cc
Dimensions:	Length: 6.7" (17cm) Height: 5.1" (13cm) Depth: 4.3" (11cm)
Input:	100-240 V 50/60 Hz 1.5 A



Alleva Medical Ltd

Suite M-Q, 12th Floor, Kings Wing Plaza 2, 1 On Kwan Street, Shatin, Hong Kong
Tel: +852 3166 7200 www.allevamedical.com

Distributed in the USA exclusively by ExtricateUSA LLC

Suite 205, 2460 General Armistead Avenue, Norristown, PA 19403
Tel: 877-312-NPWT info@extricateusa.com

All products are manufactured under US FDA, ISO 13485, and Japan PMDA registered facility of Alleva Medical Ltd.

Alleva Medical®, extriCARE®, Vide®, Maxson® are Trademarks and/or Registered Trademarks of Alleva Medical Ltd., its subsidiaries, and/or affiliates in the United States, and China and other countries. Copyright 2023.

All dimensions, weights, colours and contents may vary.



www.extricateusa.com



www.allevamedical.com

HCPCS Code: A6550



Mini Oval EC2400-MN-E Small Oval EC2400-SO-E Large Oval EC2400-LO-E Extra Large Oval EC2400-IM-E

Bandage Dressing

Features and Benefits:

- ✔ Anatomically fitted bandage dressings available in various sizes for different wound types
- ✔ Non-adhesive contact layer ensures a tight seal and allows negative pressure to reach the wound bed for optimal healing
- ✔ Patented flat, perforated drainage tube ensures even distribution of negative pressure across the wound bed
- ✔ Connects via tubing to the ExtriCARE canister for effective removal of exudate, maintaining a cleaner, healthier wound site

HCPCS Code: A6550

Foam Dressing

Features and Benefits:

- ✔ 2 step drape and easy to follow instructions
- ✔ Can customize tubing length to tailor to patient needs
- ✔ Low profile soft dome port pad
- ✔ High tensile strength foam to reduce flaking when cut to any wound size
- ✔ Large, one size fits all kit, (small and medium sizes are available if needed)
- ✔ Included ruler sticker with each kit to assist with wound Measurement



EC-Foam-L-E

Enhancing Wound Care with ExtriCARE®

The ExtriCARE® Hybrid hNPWT technology, allows the use of both traditional foam dressing kits and anatomically fitted bandage dressings with any of the ExtriCARE® NPWT pumps included the ExtriCARE® 1000, 3000, and 3600. This versatility promotes effective wound healing by providing tailored solutions for various wound types.



Bandage Dressing Specifications

Size: Extra Large, Large, Small, Mini

Size Dimensions : Mini: Pad 6 cm x 5 cm (2.36 in x 1.97 in)
Small: Pad 10.2 cm x 7.4 cm (4.02 in x 2.91 in)
Large: Pad 15.2 cm x 7 cm (5.98 in x 2.76 in)
X-Large: Pad 22 cm x 6.5 cm (8.66 in x 2.60 in)

Tubing Length: 1.5m (4.92ft)

Dressing Kit Includes:

- Anatomically fitted bandage dressing
- Connection tubing to ExtriCARE® pump
- Adhesive strips
- Ruler on packaging
- Instructions for use

Foam Dressing Specifications

Size: Large (one size fits all), Medium, Small

Size Dimensions : Small: 10cm x 7.5cm (3.94in x 2.95in)
Medium: 20cm x 13cm (7.87in x 5.11in)
Large: 25cm x 16cm (9.84in x 6.30in)

Tubing Length: 2 m (6.56 ft)

Dressing Kit Includes:

- Anatomically fitted bandage dressing
- Connection tubing to ExtriCARE® pump
- Adhesive strips
- Ruler on packaging
- Instructions for use

